



CERTIFICATE OF ANALYSIS

Name of the Product	: 100mM Acetic Acid, pH 3.0 solution	Part No.	: SS/05
Batch No.	: SS/05/0001	Batch Size	: 400 Lit.
Retest Date	: 1 Year from the date of manufacturing		
AR No.	: AR/SS/05/0001		
Packing Details	: Buffer solution packed in aseptic bag as primary packaging material and LDPE totes as secondary packaging material.		

Sr. No.	Test	Specification	Result
1.	Appearance	Clear colourless solution	Complies
2.	pH of solution	≤ 3.2	Complies
3.	Bioburden	<10 CFU/100ml	Complies
4.	Endotoxin	<1.0 EU/mL	Complies

Remark:

Storage condition of solution is preferred to store at ambient temperature.

Conclusion:

The above product complies / ~~not complies~~ as per above Specifications.



Approved By: - _____