



## CERTIFICATE OF ANALYSIS

<b>Name of the Product</b>	:	20mM Sodium phosphate, 1M sodium chloride solution	<b>Part No.</b>	:	SS/38
<b>Batch No.</b>	:	SS/38/0001	<b>Batch Size</b>	:	400 Lit.
<b>Retest Date</b>	:	1 Year from the date of manufacturing			
<b>AR No.</b>	:	AR/SS/38/0001			
<b>Packing Details</b>	:	Buffer solution packed in aseptic bag as primary packaging material and LDPE totes as secondary packaging material.			


Sr. No.	Test Parameters	Specification	Result
1.	Appearance	Clear Colourless solution	Complies
2.	pH of solution	7.2 ± 0.2	Complies
3.	Conductivity of solution	85 ± 10 mS/cm	Complies
4.	Bioburden	<10 CFU/100 mL	Complies
5.	Endotoxin	<1.0 EU/mL	Complies

**Remark:**

Storage condition of solution is preferred to store at ambient temperature.

**Conclusion:**

The above product complies / ~~not complies~~ as per above specifications.



Approved By: - \_\_\_\_\_